


DATA SHEET

**FM 3**

# MEDICAL FILTERS

**MAKE IT BECKER.**





- 1 Incl. pressure differential gauge and integrated membrane for increased protection against contamination
- 2 Incl. ball valve (97220246000)
- 3 Incl. drain flask (93130010000, 250ml)
- Dimensions in mm [inch]

	m <sup>3</sup> /h			g	mm [inch]				Order no.
	1000 mbar	-667 mbar	-917 mbar		a	b	c	d	
<b>FM 3-79 DPM</b>	15	4,5	1,5	½"	234 [9,2]	18 [0,7]	80 [3,2]	75 [3]	74007007903
<b>FM 3-83 DPM</b>	15	4,5	1,5	¾"	328 [12,9]	23 [0,9]	104 [4,1]	98 [3,9]	74007008303
<b>FM 3-93 DPM</b>	28	8,4	2,8	¾"	328 [12,9]	23 [0,9]	104 [4,1]	98 [3,9]	74007009303
<b>FM 3-113 DPM</b>	50	15	5	1"	612 [24,1]	34 [1,3]	154 [6,1]	150 [5,9]	74007011303
<b>FM 3-129 DPM</b>	75	22,5	7,5	1¼"	612 [24,1]	34 [1,3]	154 [6,1]	150 [5,9]	74007012903
<b>FM 3-133 DPM</b>	90	27	9	1½"	612 [24,1]	34 [1,3]	154 [6,1]	150 [5,9]	74007013303
<b>FM 3-143 DPM</b>	140	42	14	2"	744 [29,3]	45 [1,8]	196 [7,7]	195 [7,7]	74007014303
<b>FM 3-163 DPM</b>	170	51	17	2"	744 [29,3]	45 [1,8]	196 [7,7]	195 [7,7]	74007016303
<b>FM 3-173 DPM</b>	200	60	20	2½"	732 [28,8]	56 [2,2]	215 [8,5]	210 [8,3]	74007017303

Replacement cartridge					
for FM 3	Order no.	for FM 2 <sup>1)</sup>	Order no.	for FM <sup>1)</sup>	Order no.
<b>FM 3-79 DPM</b>	90951007903	<b>FM 2-70 DPM</b>	90951007001	<b>FM 25 (DPM)</b>	90951005000
<b>FM 3-83 DPM</b>	90951007903	<b>FM 2-80 DPM</b>	90951007001	<b>FM 50 (DPM)</b>	90951005000
<b>FM 3-93 DPM</b>	90951009303	<b>FM 2-90 DPM</b>	90951009001	<b>FM 70 (DPM)</b>	90951007000
<b>FM 3-113 DPM</b>	90951011303	<b>FM 2-110 DPM</b>	90951011001	<b>FM 90 (DPM)</b>	90951009000
<b>FM 3-129 DPM</b>	90951013303	<b>FM 2-120 DPM</b>	90951013001	<b>FM 110 (DPM)</b>	90951011000
<b>FM 3-133 DPM</b>	90951013303	<b>FM 2-130 DPM</b>	90951013001	<b>FM 120 (DPM)</b>	90951013000
<b>FM 3-143 DPM</b>	90951014303	<b>FM 2-140 DPM</b>	90951014001	<b>FM 130 (DPM)</b>	90951013000
<b>FM 3-163 DPM</b>	90951016303	<b>FM 2-160 DPM</b>	90951016001	<b>FM 140 (DPM)</b>	90951014000
<b>FM 3-173 DPM</b>	90951017303	<b>FM 2-170 DPM</b>	90951017001	<b>FM 160 (DPM)</b>	90951016000
				<b>FM 170 (DPM)</b>	90951017000

<sup>1)</sup> FM & FM 2 filter housing no longer available!  
 • Bacteria filtration rate 0.0001%, filtration efficiency (nominally) 99,9% @ 1 µm  
 • Max. operating temperature 50°C  
 • Lacquer-finished in RAL 9003 (white)

**MAKE IT BECKER.**